

Antibacterial activity

Gentamicin activity is primarily directed against aerobic, gram negative bacilli. The action against most gram positive bacteria is limited. Gentamicin is active against sensitive strains of enterococci and streptococci at concentrations which can be achieved clinically only when combined with a penicillin. Gentamicin is active *in vitro* against more than 90% of strains of *S. aureus* and 75% of *S. epidermidis*. Gentamicin has been shown to be active against most strains of the following organisms both *in vitro* and in clinical infections.

Most common susceptible pathogens

Gram positive bacteria

Staphylococcus aureus; Streptococcus pyogenes; Streptococcus pneumoniae; Streptococcus faecalis; Listeria monocytogenes

Gram negative bacteria

Citrobacter; Enterobacter; Escherichia coli Klebsiella spp.; Proteus mirabilis; Proteus vulgaris; Morganella morganii; Providencia spp.; Salmonella spp.; Serratia; Shigella spp.; Pseudomonas aeruginosa

The drug claims for gentamicin, including its mechanism of actions and antibacterial activities for the gentamicin containing bone cements have not been clinically proven. Gentamicin containing bone cements have not been shown clinically to be active against strains of the organisms indicated above.

Aminoglycosides

The aminoglycosides are a clinically important group of antibiotics that have a broad antibacterial spectrum and their action is bactericidal. The family includes streptomycin, gentamicin, tobramycin, kanamycin, amikacin and netilmicin. The aminocyclitols such as spectinomycin are closely related and have a similar mode of action. Aminoglycosides have a variety of effects within the bacterial cell but principally they inhibit protein synthesis. Another important function of the aminoglycosides is that they increase membrane leakage. Aminoglycosides have a little activity against most gram-positive bacteria, including Streptococcus pyogenes, *S. pneumoniae*, and enterococci. Although most strains of enterococci demonstrate *in vitro* resistance, some strains in this group are susceptible (*E. faecalis*). *In vitro* studies have shown that an aminoglycoside combined with an antibiotic interfering with cell-wall synthesis (e.g. penicillin, vancomycin) affects some enterococcal strains synergistically, resulting in more bactericidal effect. The β -lactam antibiotic favors the intracellular penetration of aminoglycoside. The combination of penicillin G and gentamicin results in a synergistic bactericidal effect *in vitro* against certain strains of Enterococcus faecalis. However this combination is not synergistic against other closely related organisms (e.g. *E. faecium*). Susceptibility testing and tests for antibiotic synergisms are emphasized.

Microbial resistance

Bacteria may be resistant to the antimicrobial activity of the aminoglycosides because of failure of permeation of the antibiotic, low affinity of the drug for the bacterial ribosome, or inactivation of the drug by microbial enzymes. Drug inactivation is by far the most important explanation for the acquired microbial resistance to aminoglycosides that is encountered in clinical practice. Cross-resistance between aminoglycosides may occur.

Absorption

All of the aminoglycosides are absorbed rapidly from intramuscular sites of injection. Peak concentrations in plasma occur after 30 to 90 minutes and are similar to those observed 30 minutes after completion of an intravenous infusion of an equal dose over a 30-minute period. In critically ill patients, especially those in shock, absorption of drug may be reduced from intramuscular sites because of poor perfusion. The aminoglycosides are highly polarca-tions and therefore are very poorly absorbed from the gastrointestinal tract.

Distribution

Because of their polar nature, the aminoglycosides largely are excluded from most cells, from the central nervous system, and from the eye. The apparent volume of distribution of these drugs is 25% of lean body weight and approximates the volume of extracellular fluid. Concentrations of aminoglycosides in secretions and tissues are low. High concentrations are found only in the renal cortex and in the endolymph and perilymph of the inner ear; this may contribute to the nephrotoxicity and ototoxicity caused by these drugs. Concentrations in bile approach 30% of those found in plasma as a result of active hepatic secretion, but this represents a very minor excretory route for the aminoglycosides. Penetration into respiratory secretions is poor. Diffusion into pleural and synovial fluid is relatively slow, but concentrations that approximate those in the plasma may be achieved after repeated administration. Inflammation increases the penetration of aminoglycosides into peritoneal and pericardial cavities. Concentrations of aminoglycosides in cerebrospinal fluid (CSF) that are achievable with parenteral administration of drug usually are subtherapeutic. Penetration of amino glycosides into ocular fluids is so poor that effective therapy of bacterial endophthalmitis requires periocular and intraocular injections of the drugs.

Systemic Dosing

Traditionally, the total daily dose of aminoglycosides is administered as two or three equally divided doses. Administration of the total dose once daily, however, appears to be less toxic and just as effective. Toxicity results from accumulation of drug in the inner ear and kidney. The amount of drug that accumulates increases with higher plasma concentrations and longer periods of exposure. Elimination (or washout) of aminoglycoside from these organs occurs more slowly than from plasma and is retarded by high plasma concentrations

accounting for the association between toxicity and high plasma trough concentrations. Toxicity, then, can be considered as a threshold phenomenon, more likely to occur the longer the plasma concentration exceeds a relatively safe upper limit (e.g., above a recommended trough concentration). A once daily dosing regimen, despite the higher peak concentration, provides a longer period when concentrations are below the threshold for toxicity than does a multiple-dosing regimen (12 hours versus less than 3 hours total in the example shown in the figure), accounting for its lower toxicity. Aminoglycoside bactericidal activity, on the other hand, is directly related to the concentration achieved, because aminoglycosides have concentration-dependent killing and a concentration-dependent postantibiotic effect. This enhanced activity at higher concentrations probably accounts for the equivalent efficacy of a once-daily regimen compared to a multiple-dosing regimen despite the relatively prolonged periods of time that plasma concentrations are "subtherapeutic," i.e., below the minimum inhibitory concentration (MIC).

Elimination

The aminoglycosides are excreted almost entirely by glomerular filtration, and concentrations in the urine of 50 to 200 mg/ml are achieved. A large fraction of a parenterally administered dose is excreted unchanged during the first 24 hours, with most of this appearing in the first 12 hours. The half-lives of the aminoglycosides in plasma are similar and vary between 2 and 3 hours in patients with normal renal function. Renal clearance of aminoglycosides is approximately two-thirds of the simultaneous creatinine clearance; this observation suggests some tubular reabsorption of these drugs. The concentration of aminoglycoside in plasma produced by the initial dose is dependent only on the volume of distribution of the drug. Since the elimination of aminoglycosides is almost entirely dependent on the kidney, a linear relationship exists between the concentration of creatinine in plasma and the half-life of all aminoglycosides in patients with moderately compromised renal function. In anephric patients, the half-life varies from 20 to 40 times that determined in normal individuals. Because the incidence of nephrotoxicity and ototoxicity is related to the concentration to which an aminoglycoside accumulates, it is critical to reduce the maintenance dosage of these drugs in patients with impaired renal function. Aminoglycosides are removed from the body by either hemodialysis or peritoneal dialysis.

Bibliography

Godman & Gilman's The Pharmacological Basis of Therapeutics 2005, XI Ed., Chapter 45 (Henry F. Chambers) pp.1155-1170; McGraw Hill, New York.

GENERAL PRECAUTIONS

- Prosthesis to be implanted must be compatible with the use of bone cement.
- Store in a dry place away from all sources of light at a temperature below 25°C.
- Read instruction booklet carefully.

If any form of infection should arise following surgery, patients are instructed to inform their doctors to reduce the risk of infection immediately.

Caution: NEVER add other substances or foreign bodies to UNITE® AB Bone Cement.

Caution: Bone cements reach temperatures higher than physiological temperatures during the polymerization reaction. Polymerization of the bone cement is an exothermic reaction that occurs while the bone cement is hardening in situ. The released heat may damage bone or tissue adjacent to the implant.

USER PRECAUTIONS

Contact with monomer to skin and mucous membranes should be avoided as instances of contact dermatitis have been noted in susceptible subjects. It is advised to wear a second pair of surgical gloves and carefully observe the instructions for mixing cement in order to reduce the possibility of negative reactions.

UNITE® AB Bone Cement liquid is a powerful lipid solvent. Do not come into contact with latex or rubber gloves. UNITE® AB Bone Cement should not come into contact with the gloved hand until the cement has formed into the consistency of dough (roughly 1-2 minutes after mixing). Due to the volatility and flammability of liquid monomer, it should be evaporated in a properly ventilated hood or absorbed by an inert material and transferred into a suitable container for disposal in a landfill.

The consistency of the bone cement changes in just a few minutes once the two components are mixed. Viscosity increases rapidly to form a mass which securely adheres the prosthesis to the host site.

Achievement of this state is determined by the increase in temperature of the cement. The cement cools spontaneously after a few minutes, which is the end of the reaction and time when the surgeon can release the prosthesis.

PHARMACOLOGICAL PRECAUTIONS

If gentamicin, tobramycin, or other aminoglycosides have been administered to the patient prior to surgery, monitoring of trough serum concentrations should be performed on the day before the operation. If serum concentrations exceed 1µg/ml gentamicin, tobramycin, or other aminoglycosides (amikacin, streptomycin) UNITE® AB bone cement should not be used. Simultaneous or sequential use of aminoglycosides (i.e. IV antibiotics, antibiotic beads impregnated bone cement) in patients with renal or vestibular/auditory compromise should

be avoided. Gentamicin bone cement should not be used in patients with impaired renal function (creatinine clearance less than or equal to 20 ml/min).

Monitoring

Patients receiving gentamicin bone cement should be periodically monitored with peak and trough levels of the antibiotic, serum electrolytes, serum renal function, urinalysis, and audiograms (in the elderly and/or dehydrated patient in whom there is a higher risk of adverse events associated with gentamicin use). Elderly patients may have reduced renal function that may not be evident in the results of routine screening tests, such as BUN or serum creatinine. A creatinine clearance determination may be more useful. Monitoring of renal function during treatment with aminoglycosides is particularly important in such patients. The inactivation of gentamicin and other aminoglycosides by β -lactam type antibiotics (penicillins, cephalosporins) has been demonstrated *in vitro* and in patients with severe renal impairment. Such inactivation has not been found in patients with normal renal function who have been given the drugs by separate routes of administration. Therapy with gentamicin bone cement may result in overgrowth of nonsusceptible organisms. If overgrowth of nonsusceptible organism occurs, appropriate therapy should be initiated.

Use of gentamicin bone cement should be avoided in the following situations:

- Concurrent/sequential use of other neurotoxic/ nephrotoxic antibiotics

Other antibiotics

- Cephaloridine (1st gen. β -lactams)
- Viomycin (cyclic peptides)
- Polymixins, Colistin (polypeptides)
- Cisplatin (antineoplastic agent)
- Vancomycin and Teicoplanin (glycopeptides)
- Cyclosporins (immunosuppressant)
- Amphotericin B (antifungal)

Other drugs

- Furosemide, ethacrynic acid (loop diuretics)

WARNINGS

Before, during or immediately after the use of gentamicin bone cement, consideration should be given to the administration of ototoxic or nephrotoxic drugs. This is most important in regards to elderly patients with impaired creatinine clearance and renal impairment. The use of bone cement requires a high level of coordination between the anesthetist and the surgeon. The surgeon must communicate to the anesthetist that the cement is about to be introduced during the procedure.

In some cases events defined as "bone implantation syndrome" (BCIS) may occur which are characterized by a number of clinical features that include hypoxia, hypotension, cardiac arrhythmias, increased pulmonary vascular resistance (PVR), and cardiac arrest, which must be controlled with the methods in use in modern anaesthesiology. These phenomena are commonly associated with, but is not restricted to, cemented hip arthroplasty which usually occurs at one of the five stages in the surgical procedure: femoral reaming, acetabular or femoral cement implantation, insertion of the prosthesis or joint reduction (Donaldson et al., 2009, Br J Anaesth).

Patient blood pressure should be watched during and immediately after application of bone cement. Overpressurization of the bone cement should be avoided while inserting the prostheses to minimize the possibility of pulmonary embolism.

Surgeons should be prepared and familiar with the application, properties, and handling characteristics of UNITE® AB Bone Cement to ensure a successful procedure. Curing characteristics of UNITE® AB Bone Cement may differ with their mixing technique preference and temperature, and are best determined by the surgeon. It is recommended that the surgical team perform trials prior to use in patients under the same environmental and instrumental conditions.

The operative area must be correctly ventilated as the liquid component is both flammable and volatile. The liquid monomer and vapors must never directly be exposed to flames or heated items. It has been reported that monomer vapors have ignited by use of electrocautery devices in surgical sites near newly implanted bone cement. Exercise caution while mixing the liquid and powder components of the bone cement to prevent prolonged exposure to the concentrated vapors of liquid monomer. This exposure may cause irritation of the respiratory tract, eyes, and in some cases the liver. Those with contact lenses should not prepare bone cement or be near the mixing process.

SPECIAL WARNINGS

It is not recommended to treat active infections with UNITE® AB Bone Cement.

If used under conditions not suggested it is unlikely to provide benefits to the patient and therefore may increase the risk of the development of drug-resistant bacteria.

Do not add foreign bodies or other substances (including other antibiotics) to the bone cement.

Diuretics are rarely a source of vestibulotoxicity and may be a source of hearing impairment. Diuretics may be synergistic with other aminoglycoside ototoxins such as streptomycin,

CEMENT PREPARATION PRECAUTIONS

- Do not resterilize any of the components.
- Sterility is assured only if the unit containers are not damaged or opened.
- If the powder has a yellowish or brownish color or if the liquid is syrupy, do not use product. This indicates the product has not been stored properly.
- Do not use after the expiration date since the effectiveness of the device may be compromised.
- Ensure the inner packages and components are undamaged. Powder should be consistent (no agglomerations) and not yellow or brown in color. The contents within the vial should appear as a low viscosity liquid.
- Preparation characteristics of bone cement may be affected by temperature. Temperatures of more than 23°C for the product, the preparation accessories and the environment accelerate the various stages in the preparation procedure. Lower temperatures retard the preparation stages. Prior to using UNITE® AB Bone Cement it is strongly advised to make sure that the product has been stored at a temperature of 23°C \pm 1°C for 24 hours before surgery.

CEMENT APPLICATION PRECAUTIONS

Clinical studies show the need to maintain strictly aseptic surgical procedures. Any deep infection of a surgical wound is a serious risk and will affect the successful outcome of the technique. Some infections may appear later without clinical signs years after surgery.

In order to reduce the risk of inclusion of blood and debris within the cement, and of marrow content in the vascular system, the bone cavity must be properly irrigated with Ringer or saline solutions and dried prior to the application of bone cement. It is critical to maintain the position of the prosthesis with manual pressure until the end of the polymerization process while the cement hardens; this is important to ensure ideal implantation.

SPECIAL PRECAUTIONS

Properly cemented implants are stable and last. The cement or the prosthesis (or both) may loosen or fracture due to incorrect cement insertion, trauma, or dormant infection. Follow-up on all patients regularly and in the long-term after procedure.

Extrusion of the cement beyond the area of its intended application may occur resulting in the following complications: local neuropathy; bladder fistula; hematuria; dysuria; delayed sciatic nerve entrapment from extrusion of the bone cement beyond the region of its intended application local vascular erosion and occlusion; intestinal obstruction because of the adhesion and stricture of the ileum due to heat released during the exothermic polymerization.

gentamicin, kanamycin, and neomycin. Avoid exposure to these agents in the event hearing is impaired.

The dual use of gentamicin and neuromuscular blocking agents can cause respiratory paralysis/neuromuscular blockage and may be reversed by calcium salts.

DOSAGE AND ADMINISTRATION

UNITE® AB Bone Cement (single dose) is prepared by mixing the entire contents of one packet of powder with one vial of liquid per the instructions below. Should the need arise more doses may be mixed if required.

Different lot numbers of UNITE® AB Bone Cement may be used in tandem when mixed properly.

APPLICATION INSTRUCTIONS

To improve the use of UNITE® AB Bone Cement

- Cement should be used at a temperature of 23°C \pm 1°C and relative humidity of 60%.
- Remove detritus and irrigate the bone site carefully with saline solution.
- The presence of liquid between the bone tissue and the cement should be avoided. The bone surface must be dried with gauze and/or suction catheters before and during the cementation process.
- It is important to apply an optimal thickness of bone cement. The entire stem of the prosthesis should be covered by a uniform coating of bone cement.

PREPARATION

Temperatures affect bone cement: an increase of the local environment over 23°C reduces the times, while lower temperatures increase the times shown in the table below.

Step 1

Open the package and remove internal items. Place the powder packet and vial of liquid onto a sterile working surface within the operative area.

Step 2

Break open vial and pour the liquid into a proper container for mixing. Open powder and pour over liquid. Mix the cement with a spatula from the outside of the container towards the center to minimize the presence of air bubbles. All powder must be moistened by the liquid, therefore use the spatula delicately to remove any lumps of unmoistened powder into the overall grouping of moist cement. The surgeon will determine the amount of cement needed based on the clinical application and needs.